



My Patient, Your Patient, Our Patient- Immunohematology Reference Lab

Amy Johnson, MT (ASCP)



Ultimate Goal

Carter BloodCare Mission Statement: We save lives by making transfusion possible.

This applies to every aspect from the collection of the blood, to processing, testing, distributing and everything in between.

In the Reference and Transfusion Department (R&T), one area we are specifically interested in is pre-transfusion testing. Our goal is to identify any unexpected antibodies in the patient's sample and to provide a safe blood product to the patient in a timely manner.

- ❖ What is the antibody?
- ❖ Is the antibody clinically significant?
- ❖ Can we find compatible blood?
- ❖ Does this antibody cause hemolytic disease of the newborn?
- ❖ How long is it going to take?



Let's get started

Receiving/Processing

- Lab assistants receive in specimen
- Call to notify client receipt of specimen (if reference facility)
- Centrifuge sample
- Place sample/order in designated racks



Now the fun begins

Getting started



- The next available tech will select a workup based on the order it was received and/or status or continue a workup from a previous shift
- A workup sent for antibody identification typically begins with a Direct Antiglobulin Test (DAT) and antibody screen



Testing considerations

Methodology

There is no perfect, one-size-fits-all test method.

▪ Tubes (LISS, PeG)

Pros: test different phases (IS,RT,37,IAT), resolve discrepancies/rouleaux, specialized testing
Cons: tech/technique dependent, labor intensive, subjective, unstable (must be read immediately)

▪ Gel

Pros: small sample/reagent volume, less subjective interpretation, stable, can be automated
Cons: time limitation, sensitivity (not enough/too much), equipment requirements/maintenance

▪ Solid Phase (ECHO)

Pros: automated, stable, standardization, detect Kidd antibodies well, Knops system not detected
Cons: larger volume requirement, unable to choose selected cells, sensitivity, time limitation



Collaboration

▪ General Considerations

- Methodology used at facility (gel, tubes, solid phase)
- Results from testing performed at facility
 - +DAT?
 - Anagram results? All cells pos? 1 cell? Incompat xm?
- Previous history check
 - WAA/CAA
 - Allo Ab? Ex: -K, -Js(b)
- Additional Details
 - Recently transfused? BMT?
 - Diagnosis?
 - Multiple Myeloma, ITP, Leukemia, Sickle Cell
 - Ethnicity?
 - Medication?
 - WinRho, Daratumumab (DARA), anti-microbial (Cephalosporin)



Conclusion

	Rh-Hr			Kell			Duffy		Kidd		Lewis		P		MN		Low Ph EL	LW	2X Peg Ads plasma *			
	D	C	c	E	e	f	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M				N	S	s
1	+	+	0	0	+	0	+	+	+	0	+	0	+	0	+	0	0	+	+	2+	0✓	0✓
2	+	0	+	+	0	0	0	+	+	+	0	+	0	+	+	0	0	+	+	2+	0✓	0✓
3	0	0	+	+	+	+	+	+	0	+	0	+	0	0	+	+	0	+	+	2+	0✓	0✓
*			%			%	0				0									2+		0✓

*W035215235237-E (adsorbing cell)



Warm auto Ab with no underlying Allo antibodies



Case study 1

- 37 year old pregnant, Hispanic female
- Presented to hospital for delivery with a 3g/dl Hgb
- Referred to our facility by a reference account

	Rh-Hr			Kell			Duffy		Kidd		Lewis		P		MN		LSS 37*	LSS IAT	Peg IAT	Gel		
	D	C	c	E	e	f	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M					N	S
1	+	+	0	0	+	0	+	+	+	0	+	0	+	0	+	0	0	+	+	1+	2+s	2+
2	+	0	+	+	0	0	0	+	+	0	+	0	+	0	+	0	0	+	+	1+s	3+	2+
3	0	0	+	+	+	+	+	+	0	+	0	+	0	0	+	+	0	0	0	1+	2+s	3+
AC																		0	0✓	0✓	0	0



Next step

Name: _____		Patient ID: _____		Code: _____		Room: _____									
Date: _____		Order: _____		Order: _____		Order: _____									
ABO/TYPING				ADDITIONAL TYPING				GAT							
A	B	AB	O	A	B	AB	O	Anti-A	Anti-B	Anti-AB	Anti-O	Anti-A	Anti-B	Anti-AB	Anti-O
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
ANTIBODY SCREEN				ADDITIONAL RED BLOOD CELL PHENOTYPES				Cell Treatment							
Anti-A	Anti-B	Anti-AB	Anti-O	Anti-A	Anti-B	Anti-AB	Anti-O	Anti-A	Anti-B	Anti-AB	Anti-O	Anti-A	Anti-B	Anti-AB	Anti-O
+	+	+	+	+	+	+	+	+	+	+	+	+	+	+	+
MISCELLANEOUS TESTING															
RESULTS / CONCLUSIONS:															



Enzymes/Adsorption

	Rh-Hr			Kell			Duffy		Kidd		Lewis		P		MN		PeG IAT	JM DTT	Ficoll 37/IAT	#	
	D	C	c	E	e	f	K	k	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M					N
1	+	+	0	0	+	0	0	+	0	+	0	+	0	0	+	0	0	+		0	✓
2	+	0	+	0	0	0	+	0	0	+	0	+	0	+	0	+	0	+		0	✓
3	0	0	+	0	+	+	+	0	+	0	+	0	0	+	+	+	0	+		0	✓
*W03521 6084830	0					0			0								0	3+	3+	0/0	0
W03521 5246554	(Dib=)																	2+s			0
O cord																		3+			0

*phenotypically matched
2x PeG adsorbed plasma w/ *



The Blood Group Antigen Facts Book

Ficoll/ Papain	Trypsin	α-Chymotrypsin	200 mM DTT/AET	Possible specificity
Negative	Negative	Negative	Positive	Bp ^a , Ch ^a , Ig, JG5
Negative	Negative	Negative	Negative	PN, BMH ^a
Negative	Negative	Positive	Positive	HL, PL, E ^a , F ^a , Ge ^a , Ge ^b , Ge ^c
Negative	Positive	Negative	Positive	Y ^a , Fy ^a , Fy ^b
Variable	Positive	Negative	Positive	S, s
Variable	Positive	Negative	Weak or negative	Y ^a
Negative	Positive	Positive	Positive	Ro ^a , S ^a
Positive	Negative	Negative	Weak or negative	LU, MEX ^a
Positive - Papain Weak or negative - Ficoll	Negative	Negative	Negative	RN ^a
Positive	Negative	Weak	Negative	DOG
Positive	Positive	Negative	Weak	CRD ^a
Positive	Positive	Negative	Positive	Some De (L ^a) range
Positive	Positive/weak	Positive	Negative	SW ^a
Positive	Positive/weak	Positive	Positive	SC
Positive	Positive ^a	Positive ^a	Negative	KEL ^a (except KAL ^a , which is transfer-inhibitable)
Positive	Positive	Positive	Positive	ABO, E ^a , R ^a , LL, PL ^a , R ^b , S ^a , F ^a , PL, H ^a , H ^b , K ^a , K ^b , P ^a , Ge ^a , Ge ^b , Ge ^c , Lu, P ^a , F ^a , F ^b , F ^c , S ^a , S ^b , S ^c , S ^d , S ^e , S ^f , S ^g , S ^h , S ⁱ , S ^j , S ^k , S ^l , S ^m , S ⁿ , S ^o , S ^p , S ^q , S ^r , S ^s , S ^t , S ^u , S ^v , S ^w , S ^x , S ^y , S ^z , S ¹ , S ² , S ³ , S ⁴ , S ⁵ , S ⁶ , S ⁷ , S ⁸ , S ⁹ , S ¹⁰ , S ¹¹ , S ¹² , S ¹³ , S ¹⁴ , S ¹⁵ , S ¹⁶ , S ¹⁷ , S ¹⁸ , S ¹⁹ , S ²⁰ , S ²¹ , S ²² , S ²³ , S ²⁴ , S ²⁵ , S ²⁶ , S ²⁷ , S ²⁸ , S ²⁹ , S ³⁰ , S ³¹ , S ³² , S ³³ , S ³⁴ , S ³⁵ , S ³⁶ , S ³⁷ , S ³⁸ , S ³⁹ , S ⁴⁰ , S ⁴¹ , S ⁴² , S ⁴³ , S ⁴⁴ , S ⁴⁵ , S ⁴⁶ , S ⁴⁷ , S ⁴⁸ , S ⁴⁹ , S ⁵⁰ , S ⁵¹ , S ⁵² , S ⁵³ , S ⁵⁴ , S ⁵⁵ , S ⁵⁶ , S ⁵⁷ , S ⁵⁸ , S ⁵⁹ , S ⁶⁰ , S ⁶¹ , S ⁶² , S ⁶³ , S ⁶⁴ , S ⁶⁵ , S ⁶⁶ , S ⁶⁷ , S ⁶⁸ , S ⁶⁹ , S ⁷⁰ , S ⁷¹ , S ⁷² , S ⁷³ , S ⁷⁴ , S ⁷⁵ , S ⁷⁶ , S ⁷⁷ , S ⁷⁸ , S ⁷⁹ , S ⁸⁰ , S ⁸¹ , S ⁸² , S ⁸³ , S ⁸⁴ , S ⁸⁵ , S ⁸⁶ , S ⁸⁷ , S ⁸⁸ , S ⁸⁹ , S ⁹⁰ , S ⁹¹ , S ⁹² , S ⁹³ , S ⁹⁴ , S ⁹⁵ , S ⁹⁶ , S ⁹⁷ , S ⁹⁸ , S ⁹⁹ , S ¹⁰⁰
Positive	Positive	Positive	Enhanced	S ^a

*Kell blood group system antigens are sensitive to treatment with a mixture of trypsin and α-chymotrypsin. DTT may be variable.



Confirmation



➤ Type patient for Ge2 Ag



➤ Test Ge:-2 cells against patient's plasma



Blood selection

GE | Gerbich Blood Group System 511

Effect of enzymes and chemicals on Ge2 antigen on intact RBCs

Ficin/Papain	Sensitive
Trypsin	Sensitive
α -Chymotrypsin	Weakened
Sialidase	Variable
DTT 200mM	Variable (thus variable to WARM™ and ZZAP)
Acid	Resistant

In vitro characteristics of alloanti-Ge2

Immunoglobulin class	Dually IgG
Optimal technique	IAT
Complement binding	Yes; some hemolytic

Clinical significance of alloanti-Ge2

Transfusion reaction	No to moderate/immmediate/delayed
HDPN	Positive DEX, but no clinical HDPN

Autoanti-Ge2

Yes; detects a determinant on GPC.

Comments

Alloanti-Ge2 can be made by individuals with Yes, Gerbich or Leach phenotypes; detects an antigen on GPC, and may be naturally-occurring. The reciprocal gene to GPC, *hcr* encodes two copies of amino acids encoded by *hcr* 2.

Siblings of patients with anti-Ge2 should be tested for compatibility, and the patient urged to donate blood for cryogenic storage when in her clinical state permits.

Read HE and Larissa/Farrell, C. The Blood Group Antigen Facts Book. San Diego, CA: American Press, 3rd Edition, 2012.



- Import/ARDP
- Type siblings
- Autologous donation when available



22

What's new and exciting?

- **Molecular**
 - Screening donors
 - Testing patients
- **Daratumumab (DARA)**
 - Monoclonal anti-CD38
 - Promising for Multiple Myeloma patients



23

It takes a TEAM effort

- **Everyone plays an integral role in getting the proper blood to the patient**
 - Phlebotomist
 - Referring facility
 - Physician
 - Nurse
 - Lab
 - Family



24

Questions?



Photo credit: Bernard Simon



References

1. Reid ME and Lomas-Francis C. The Blood Group Antigen Facts Book. San Diego, CA: American Press, 3rd Edition, 2012


