



To: Carter BloodCare Customers
From: Hospital Relations
Date: October 22, 2020
Re: Imported Convalescent Plasma Product Codes

Carter BloodCare has imported COVID-19 convalescent plasma (CCP) units to support the increase in community usage. These units will be distributed immediately to fill orders as needed.


The imported CCP units are labeled as plasma frozen within 24 hours (PF24) with an attached orange manila tag identifying the unit(s) as convalescent plasma. This type of labeling is acceptable per ISBT 128 and FDA guidelines. Examples of the product codes and tag are provided for your reference.

Impact to Client

- Imported PF24 labeled CCP units will automatically have an associated plasma fee on the pack list; however, a credit memo will be issued immediately.
- Notify individuals responsible for billing to apply a \$0.00 charge for PF24 labeled CCP units as CCP does not have an associated fee.
- Imported apheresis CCP collections will be classified as Part A0, B0, C0 or D0.
- Imported convalescent plasma units could be labeled as high or low titer. This titer test does not meet the December 1st FDA titer testing requirements.
- Continue to refer recovered COVID-19 patients to carterbloodcare.org to consider donating convalescent plasma.


Thank you for your dedication and efforts to help meet patient needs in our community.

Identifying Tag



ISBT Product Codes

VOLUNTEER DONOR



E1624VA0


**APHERESIS
PLASMA
DIVIDED
FROZEN WITHIN 24 HOURS
AFTER PHLEBOTOMY**

208 mL containing approx 28 mL ACD – A
Store at – 18 C or colder

Part A0

Rx Only

VOLUNTEER DONOR




E1624VB0

**APHERESIS
PLASMA
DIVIDED
FROZEN WITHIN 24 HOURS
AFTER PHLEBOTOMY**

207 mL containing approx 31 mL ACD – A
Store at – 18 C or colder

Part B0

VOLUNTEER DONOR



E1624VC0


**APHERESIS
PLASMA
DIVIDED
FROZEN WITHIN 24 HOURS
AFTER PHLEBOTOMY**

208 mL containing approx 30 mL ACD – A
Store at – 18 C or colder

Part C0

Rx Only

VOLUNTEER DONOR



E1624VD0

**APHERESIS
PLASMA
DIVIDED
FROZEN WITHIN 24 HOURS
AFTER PHLEBOTOMY**

195 mL containing approx 26 mL ACD – A
Store at – 18 C or colder

Part D0